



**Fraunhofer** Institut  
Angewandte  
Polymerforschung

Indication of loading /  
coupling of colloidal  
carriers

Synthetic polymer colloids are applied as carriers of biomolecules. Their surface properties often change dramatically by coupling enzymes or by adsorption of proteins. That is reflected by comparisons of electrophoretic mobility at variable pH.

**Fraunhofer-Institut für  
Angewandte Polymerforschung**

Wissenschaftspark Golm  
Geiselbergstraße 69  
14476 Potsdam  
Germany

Phone +49(0)331/568-10  
Fax +49(0)331/568-3000  
E-mail [info@iap.fraunhofer.de](mailto:info@iap.fraunhofer.de)  
[www.iap.fraunhofer.de](http://www.iap.fraunhofer.de)